

RMRC Genotyping Protocol-PCR

No.	RMRC13404	Strain Name	C57BL/6- <i>Tpt1</i> ^{em1Yen} /Narl
-----	-----------	-------------	--

A. Primer

	Primer name	Sequence
1	Tpt1 e2 GT Fw	5'-tctatatgaggttggggagcgc-3'
2	B6 tpt1 E12 Rv	5'-gatctttagatgtcggagaacagct -3'
3	Tpt1 V12 Rv	5'-gatctttagatgtcggagaacagca -3'
4	Internal control Fw	5'-gccacttgctaaagcacctc -3'
5	Internal control Rv	5'-tgcacatggctctcttctcc -3'

B. Genomic extraction :

- ☐ Phenol-Chloroform extraction
 ☐ commercial extraction kit _____

☒ Lysis buffer _____
 ☐ other _____

C. Reaction Conditions: for Tpt1 wt allele

Component	Concentration		Volume		Step	Tempera- ture, °C	Time, min	Number of Cycles
H ₂ O			1	μl	1	94	3	1
Ampliqon Taq DNA Polymerase 2x Master Mix RED	2X		12.5	μl	2	94 68 (-0.5°C per cycle) 72	0.5 0.5 0.5	30
Primer 1	Tpt1 e2 GT Fw	10 μM	0.375	μl	3	94	20sec	23
Primer 2	B6 tpt1 E12 Rv	10 μM	0.375	μl		62	0.5	
Primer 4	Internal control Fw	10 μM	0.375	μl		72	35sec	
Primer 5	Internal control Rv	10 μM	0.375	μl				
5X CES	5X		5	μl	4	72	5	1
DNA	--		5	μl				
Total			20	μl	5	4	pause	

C. Reaction Conditions: for Tpt1 E12V allele

Component	Concentration		Volume		Step	Tempera- ture, °C	Time, min	Number of Cycles
H ₂ O			1	μl	1	94	3	1
Ampliqon Taq DNA Polymerase 2x Master Mix RED	2X		12.5	μl	2	94 68 (-0.5°C per cycle) 72	0.5 0.5 0.5	30
Primer 1	Tpt1 e2 GT Fw	10 μM	0.375	μl	3	94	20sec	23
Primer 3	tpt1 V12 Rv	10 μM	0.375	μl		62	0.5	
Primer 4	Internal control Fw	10 μM	0.375	μl		72	35sec	
Primer 5	Internal control Rv	10 μM	0.375	μl				
5X CES	5X		5	μl	4	72	5	1
DNA	--		5	μl				

The information of Taq, including company name and catalog number :

Taq DNA Polymerase 2x Master Mix RED, 2 mM MgCl₂; Ampliqon #A190303

D. PCR Product

	Primer set		Product size	Target region
1	Tpt1 e2 GT Fw	B6 tpt1 E12 Rv	350 bp	
2	Tpt1 e2 GT Fw	tpt1 V12 Rv	350 bp	
3	Internal control Fw	Internal control Rv	492 bp	

E. Gel Photo

